


CM1200A

12-channel Electrocardiograph

Features

- 8.4 inch LED touch screen
- Multi-language operation and information input
- 12-lead acquisition and amplification, 12 ECG waveforms recording simultaneously
- 1:1 ratio on-screen diagnose
- Manual/ Auto/ Rhythm working modes selectable
- AC/DC power supply, built-in battery, support 4 hours working time
- 300 ECG records storage and review, built-in SD card(optional) can save 40000 cases
- Support 210mm, 216mm roll and Z-folded paper
- Net-grid on screen for easier observation

Interface



Various printing modes



Printing mode: 12X1, 6X2+1R, 3X4+1R

COMEN



CM1200A

12-channel Electrocardiograph



Shenzhen Comen Medical Instruments Co., Ltd.

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district, Guangming District, Shenzhen, 518106, Guangdong, China.
Tel:+86-755-2640 8879 400-700-9488 Fax:+86-755-2643 1232
Website: en.comen.com E-mail: info@szcomen.com

P/N: EN-CM1200A-4P-20200302-V1.3





CM1200A

12-channel Electrocardiograph



8.4 inch LED touch screen



Alphanumeric keyboard.




Hidden handle



Memory 300 ECG records;
USB for data transmission;
Support PC-ECG software;
Support external printer connection



Clear professional report



12X1 + average template



Supports USB printing, A4 paper printing without computer, long-term preservation and saving cost.
Supports USB flash disk back-up and recovery of patient data, internal storage expanded; Supports connction to computer for preserving and reading patient data.
Supports preservation and copy of cases in picture format.