

AN ISO CERTIFIED COMPANY



SEMS SURGICAL INDUSTRIES

MANUFACTURERS & IMPORTERS OF HOSPITAL FURNITURE,
NICU, ICU, OT EQUIPMENTS & ETP



SEMS SURGICAL INDUSTRIES

HOSPITAL BEDS

SEMS S101: MOTORIZED REMOTE OPERATED ICU BED

- FOUR SECTIONS MS MATTRESS PLATFORM.
- MOTORIZED BACKREST, KNEE REST, TB/RTB & HEIGHT ADJUSTMENT
- ABS MOULDED SPLIT TYPE SAFETY SIDE RAILS (SWING TYPE)
- ABS MOULDED HEAD & FOOT PANELS, EASILY REMOVABLE WITHOUT USING ANY TOOLS.
- 125 MM DIAMETER SINGLE WHEEL HIGH GRADE SYNTHETIC BODY CASTORS WITH DIAGONAL LOCKING MECHANISM

OPTIONAL-

- SS HEAVY DUTY I.V. ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.

SIZE	80" (L) X 36" (W)
HEIGHT RANGE	20" - 36"
FRAMEWORK MATERIAL & FINISH	MILD STEEL PRE-TREATED & EPOXY POWDER COATED



SEMS S102: FOUR FUNCTION ICU BED WITH ABS PANEL

- FOUR SECTIONS MS MATTRESS PLATFORM.
- MANUAL KEY OPERATED BACKREST, KNEE REST, TB/RTB & HEIGHT ADJUSTMENTS.
- ABS MOULDED SPLIT TYPE SAFETY SIDE RAILS (SWING TYPE)
- ABS MOULDED HEAD & FOOT PANELS, EASILY REMOVABLE WITHOUT USING ANY TOOLS.
- 125 MM DIAMETER SINGLE WHEEL HIGH GRADE SYNTHETIC BODY CASTORS WITH DIAGONAL LOCKING MECHANISM

OPTIONAL :-

- S.S DETACHABLE COLLAPSIBLE SIDE-RAILS.
- S.S. HEAD AND FOOT BOWS.
- S.S. HEAVY DUTY I.V. ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.



SIZE	80" (L) X 36" (W)
HEIGHT RANGE	20" - 36"
FRAMEWORK MATERIAL & FINISH	MILD STEEL PRE-TREATED & EPOXY POWDER COATED

SEMS S103: THREE FUNCTION ICU BED

- FOUR SECTIONS MS MATTRESS PLATFORM.
- MANUAL KEY OPERATED BACKREST, KNEE REST, TB/RTB & HEIGHT ADJUSTMENTS.
- ABS MOULDED SPLIT TYPE SAFETY SIDE RAILS (SWING TYPE)
- ABS MOULDED HEAD & FOOT PANELS, EASILY REMOVABLE WITHOUT USING ANY TOOLS.
- 125 MM DIAMETER SINGLE WHEEL HIGH GRADE SYNTHETIC BODY CASTORS WITH DIAGONAL LOCKING MECHANISM

OPTIONAL :-

- S.S DETACHABLE COLLAPSIBLE SIDE-RAILS.
- S.S. HEAD AND FOOT BOWS.
- S.S. HEAVY DUTY IV ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.

SIZE	80" (L) X 36" (W) X 24" (H)
FRAMEWORK MATERIAL & FINISH	MILD STEEL PRE-TREATED & EPOXY POWDER COATED



SEMS S104: FOWLER BED (DELUXE)

- MANUAL KEY OPERATED BACKREST & KNEE REST ADJUSTMENT1
- ABS MOULDED HEAD & FOOT PANELS, EASILY REMOVABLE WITHOUT USING ANY TOOLS
- 100 MM DIAMETER SINGLE WHEEL HIGH GRADE SYNTHETIC BODY CASTORS WITH DIAGONAL LOCKING MECHANISM.

OPTIONAL :-

- SS HEAVY DUTY IV ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.
- S.S DETACHABLE COLLAPSIBLE SIDE-RAILS.



Size	80" (L) X 36" (W) X 24" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated



SEMS SURGICAL INDUSTRIES

HOSPITAL BEDS

SEMS S105: FOWLER BED

- MANUAL KEY OPERATED BACKREST & KNEE REST ADJUSTMENT
- S.S./M.S. HEAD & FOOT PANELS, EASILY REMOVABLE.

OPTIONAL :-

- SS HEAVY DUTY IV ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.
- S.S DETACHABLE COLLAPSIBLE SIDE-RAILS.

Size	80" (L) X 36" (W) X 24" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated



SEMS S106: SEMI FOWLER BED (DELUXE)

- MANUAL KEY OPERATED BACKREST ADJUSTMENT
- ABS MOULDED HEAD & FOOT PANELS, EASILY REMOVABLE WITHOUT USING ANY TOOLS
- 100 MM DIAMETER SINGLE WHEEL HIGH GRADE SYNTHETIC BODY CASTORS WITH DIAGONAL LOCKING MECHANISM.
- SINGLE SECTIONS MS MATTRESS PLATFORM.

OPTIONAL:

- S.S DETACHABLE COLLAPSIBLE SIDE-RAILS.
- SS HEAVY DUTY IV ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.

Size	77" (L) x 36" (W) x 24" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated

SEMS S107: SEMI FOWLER BED

- MANUAL KEY OPERATED BACKREST ADJUSTMENT
- S.S./M.S. HEAD & FOOT PANELS, EASILY REMOVABLE
- SINGLE SECTIONS MS MATTRESS PLATFORM.

OPTIONAL:

- S.S DETACHABLE COLLAPSIBLE SIDE-RAILS.
- S.S. HEAVY DUTY I.V. ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.

Size	77" (L) x 36" (W) x 24" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated



SEMS S108: PLAIN BED (PREMIUM)

- ABS MOULDED HEAD & FOOT PANELS, EASILY REMOVABLE WITHOUT USING ANY TOOLS
- URINE BAG HOLDER
- SINGLE SECTIONS MS MATTRESS PLATFORM.

OPTIONAL:

- S.S DETACHABLE COLLAPSIBLE SIDE-RAILS.
- S.S. HEAVY DUTY I.V. ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.

Size	72" (L) x 36" (W) x 24" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated

**SEMS S109: PLAIN BED**

- S.S./M.S. HEAD & FOOT PANELS, EASILY REMOVABLE.
- URINE BAG HOLDER
- SINGLE SECTIONS MS MATTRESS PLATFORM.

OPTIONAL:

- S.S DETACHABLE COLLAPSIBLE SIDE-RAILS.
- S.S. HEAVY DUTY I.V. ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.

Size	72" (L) x 36" (W) x 24" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated

**SEMS S110: PEDIATRIC BED DELUX**

- MANUAL KEY OPERATED BACKREST ADJUSTMENT
- S.S. TUBULAR, TELESCOPIC, COLLAPSIBLE SIDE RAILS
- S.S./ M.S. HEAD & FOOT PANELS
- URINE BAG HOLDER

OPTIONAL :-

- S.S. HEAVY DUTY IV ROD (TWO HOOKS) FOUR SIDE FIXING FACILITY.

Size	60" (L) x 30" (W) x 24" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated

**SEMS S111: PEDIATRIC BED**

- CRCA SHEET TOP WITH PROPER SUPPORT.
- HEAD AND LEG BOWS OF EQUAL HEIGHT WITH TWO VERTICAL AND ONE HORIZONTAL SUPPORT.
- LEG FITTED RUBBER SHOES.
- TWO SIDE FULL LENGTH DETACHABLE RAILS..

Size	60" (L) x 30" (W) x 24" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated

**SEMS S112: RECLINING CHAIR CUM ATTENDANT BED**

- CRCA TOP WITH TUBULAR FRAMEWORK.
- FIXED UPHOLSTERED TOP.

Size	30" - 72" (L) x 24" (W) x 18" (H)
Framework Material & finish	Mild Steel Pre-Treated & Epoxy Powder Coated

**SEMS S113: ATTENDANT BED**

- S.S./M.S. HEAD & FOOT PANELS.
- SINGLE SECTIONS MS MATTRESS PLATFORM.



SIZE	72" (L) X 24" (W) X 18" (H)
FRAMEWORK MATERIAL & FINISH	MILD STEEL PRE-TREATED & EPOXY POWDER COATED

**SEMS S114: LABOUR BED**

- EASILY REMOVABLE ABS MOULDED HEAD AND FOOT PANELS.
 - MINIMUM HEIGHT OF 530MM FOR THE CONVENIENCE OF PATIENTS TO GET IN AND OUT OF THE BED.
 - 125MM DIA. CENTRAL LOCKING CASTERS FOR BETTER STABILITY.
 - REMOVABLE MIDDLE AND LEG MATTRESS PLATFORM SECTION GIVES CLEAR ACCESS TO BED FROM TOP AND BOTTOM FOR EFFECTIVE CLEANING.
 - ANGULAR AND HEIGHT ADJUSTABLE.
 - HEIGHT ADJUSTABLE LITHOTOMY CRUTCHES.
- OPTIONAL :-**
- S.S. HEAVY DUTY IV ROD.
 - ABS SIDE RAILING.

SEMS S115: LABOUR TABLE

- S.S./M.S. RECTANGULAR TUBE TOP FRAME.
- S.S./M.S. BOTTOM FRAME.
- BACKREST ON RATCHET MECHANISM.
- TRENDLENBURG/REVERSE TRENDLENBURG ON SCREW MECHANISM.
- MIDDLE SECTION HAVING "U" NOTCH, WITH BOWL.
- SLIDING FOOT END.
- 35MM THICK MATTRESS.

OPTIONAL :-

- STEP STOOL (CODE: 603).
- LITHOTOMY CRUTCHES
- S.S./M.S. RAILING ON BOTH SIDES.





SEMS SURGICAL INDUSTRIES

HOSPITAL COUCH

SEMS S201: EXAMINATION COUCH

- CRCA FRAME WITH CABINETS.
- FIXED UPHOLSTERED TOP IN TWO SECTIONS.
- HEADREST ADJUSTABLE ON GAS SPRING.
- THREE CABINETS AND THREE DRAWERS.
- INBUILT SLIDING STEP STOOL
- B.P. APPARATUS TRAY.



SIZE	72" (L) X 20" (W) X 32" (H)
FRAMEWORK MATERIAL & FINISH	MILD STEEL PRE-TREATED & EPOXY POWDER COATED

SEMS S202: EXAMINATION TABLE



- CRCA TUBULAR FRAMEWORK.
- ONE DRAWER AND ONE CABINET.

SIZE	72" (L) X 20" (W) X 30" (H)
FRAMEWORK MATERIAL & FINISH	MILD STEEL PRE-TREATED & EPOXY POWDER COATED



SEMS S203:
Obstetrics-Gynaec table



SEMS S204:
Examination couch "U" cut



SEMS S205:
Examination table

HOSPITAL STRETCHER

SEMS S301: EMERGENCY & RECOVERY TROLLEY



- TWO SECTION TOP WITH REMOVABLE X-RAY PERMEABLE STRETCHER
- HEIGHT ADJUSTMENT BY SMOOTH CRANK MECHANISM
- RAISED BACKREST ON RATCHET
- HIGH GRADE SYNTHETIC BODY CASTORS WHEELS (125 MM) TWO WITH BRAKES & TWO WITHOUT BRAKES
- PUSH-PULL HANDLE & SAFETY BELTS
- PRETREATED AND EPOXY POWDER COATED
- SWING AWAY TYPE SIDE RAILS
- SUPPLIED IN SEMI KNOCK DOWN CONDITION

OPTIONAL:-

STAINLESS STEEL IV POLE ALONG WITH PROVISION AT 4 POSITIONS



SEMS SURGICAL INDUSTRIES

HOSPITAL STRETCHER



SEMS S302
FOLDING STRETCHER



SEMS S303
STRETCHER ON TROLLEY



SEMS S304
STRETCHER ON TROLLEY
(BACKREST ADJUSTABLE)

TROLLEYS



SEMS S401
ECG TROLLEY



SEMS S402
INSTRUMENT TROLLEY



SEMS S403
Defibrillator TROLLEY



SEMS S404
MEDICINE CRASH CART
CUM DRUG TROLLEY



SEMS S405
MEDICINE CRASH CART



SEMS S406
MONITOR TROLLEY



SEMS S407
DRESSING TROLLEY



SEMS S408
MAYO TROLLEY



SEMS S409
C-CURVE TROLLEY



SEMS SURGICAL INDUSTRIES

TROLLEYS'



SEMS S410
SOLID LENIN TROLLEY



SEMS S411
LENIN TROLLEY



SEMS S412
DRUM TROLLEY



SEMS S413
CYLINDER TROLLEY

BED-SIDE LOCKERS



SEMS S501
NORMAL LOCKER



SEMS S502
3-TOP LOCKER



SEMS S503
DELUX LOCKER (S.S. TOP)



SEMS S504
DELUX LOCKER (ABS TOP)



SEMS S505
PREMIUM LOCKER



SEMS S601
REVOLVING STOOL



SEMS S602
ATTENDANT STOOL



SINGLE STEP



DOUBLE STEP

SEMS S603 FOOT STEP

STOOLS

OVER BED TABLES



SEMS S701
**OVER BED TABLE
(FIXED TOP)**



SEMS S702
**OVER BED TABLE
(ADJUSTABLE TOP)**



SEMS SURGICAL INDUSTRIES

OTHER EQUIPMENTS



SEMS S801
BED SIDE SCREEN



SEMS S802
BOWL STAND



SEMS S803
KICK BUCKET



SEMS S804
CRIB CARRIER



SEMS S805
SCRUBBER



SEMS S806
WHEEL CHAIR



SEMS S807
ALMIRAH



SEMS S808
STAFF LOCKER



SEMS S809

SALINE STAND



SEMS S810

BOOK SHELF





SEMS SURGICAL INDUSTRIES

OTHER EQUIPMENTS



SEMS S811

BLOOD DONATION CHAIR



SEMS S812

DOG TABLE



SEMS S813

MONITOR STAND



SEMS S814

ELECTRIC EXAMINATION TABLE



SEMS S815

LOBBY BENCH



SEMS S816

TILT TABLE



SEMS S817

SUSPENSION BED



SEMS S818

CANTEEN TABLE



SEMS S819

FOOD DISTRIBUTION TROLLEY



SEMS SURGICAL INDUSTRIES

ACCESSORIES




SEMS S901

MATTRESS FOR PLAIN BED/ATTENDANT BED

SEMS S902

**TWO SECTION MATTRESS FOR
SEMI-FOWLER BED**



SEMS S903

**FOUR SECTION MATTRESS FOR
FOWLER BED & ICU BED**

SEMS S904

**MATTRESS FOR
LABOR BED/LABOR TABLE**



SEMS S905

**HEAD SPLIT MATTRESS FOR
EXAMINATION COUCH &
EXAMINATION TABLE**



SEMS S906

**TWO SECTION , U-CUT MATTRESS
FOR GYNAEC TABLE AND
EXAMINATION COUCH**




SEMS SURGICAL INDUSTRIES

ACCESSORIES



SEMS S907 ABS PANELS'



SEMS S908 ABS SIDE RAILINGS



SEMS S909 S.S./M.S. SIDE RAILINGS



**SEMS S910
CASTER WHEEL**



**SEMS S913
OVER-BED TABLE TOP**



**SEMS S911
IV RODS
(TWO HOOK/FOUR HOOK)**



**SEMS S912
S.S./M.S. PANEL**



**SEMS S914
ACTUATORS' & HYDRAULICS**



SEMS S914





SEMS SURGICAL INDUSTRIES

O.T. EQUIPMENT



OT-101
C-ARM COMPATIBLE OT TABLE



OT-102
HYDRAULIC OT TABLE



OT-103

DOUBLE DOME CEILING LED LIGHT 4 + 4



Intensity	1,20,000+1,20,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	4300 - 5200 K
Colour Rendering Index	95 Ra
Focusing	Adjustable
LED Average Life	>50000 Hrs.
Number of LED	76+76
Diameter of Light	650 mm
Intensity Control	10-99 Digital Touch
Power Supply	220V/50Hz AC drive by DC 24 volt
Power Consumption	80+80 Watts
Special Function	remote operated

Intensity	1,50,000+1,50,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	4300 - 5200 K
Colour Rendering Index	95 Ra
Focusing	Fix
LED Average Life	>50000 Hrs.
Number of LED	54+54
Diameter of Light	550 mm
Intensity Control	10-99 Digital Touch
Power Supply	220V/50Hz AC drive by DC 24 volt
Power Consumption	60 Watts
Special Function	remote operated

OT-104

DOUBLE DOME CEILING LED LIGHT



Intensity	1,50,000 Lux
Spot Diameter	100 - 150 mm
Depth Illumination	10 - 15 cm
Colour Temp.	4300 - 5200 K
Colour Rendering Index	95 Ra
Focusing	Adjustable
LED Average Life	>50000 Hrs.
Number of LED	95
Diameter of Light	650 mm
Intensity Control	10-99 Digital Touch
Power Supply	220V/50Hz AC drive by DC 24 volt
Power Consumption	100 Watts
Special Function	remote operated

OT-105

SINGLE DOME CEILING LED LIGHT



Intensity	60,000	60,000
Spot Diameter	100-120mm	100-120mm
Depth Illumination	10 - 12 mm	10 - 12 mm
Colour Temp.	3800 - 4500 K	3800 - 4500 K
Colour Rendering Index	91 Ra	91 Ra
Focusing	FIX	FIX
LED Average Life	>50000 Hrs.	>50000 Hrs.
Number of LED	19	21
Diameter of Light	160 mm	160 mm
Intensity Control	1-9 Digital Touch	1-9 Digital Touch
Power Supply	220V/50Hz AC drive by DC 24 volt	220V/50Hz AC drive by DC 24 volt
Power Consumption	20 Watts	20 Watts
LED Colour	Yellow & White Combination	Yellow & White Combination



OT-106

EXAMINATION LIGHT



SEMS SURGICAL INDUSTRIES

AWARDS'

EXCELLENCE IN HOSPITAL EQUIPMENT & FURNITURE BY MY FM HEALTHCARE EXCELLENCE AWARDS



GAURAVSHALI MADHYA PRADESH AWARD FOR NON STOP WORK & CONTRIBUTION DURING COVID PANDEMIC

INDUSTRY LEADERS AWARDS 2022 FOR MOST TRUSTED HOSPITAL FURNITURE MANUFACTURER IN CENTRAL INDIA BY BRAND EMPOWER (P) LIMITED



CERTIFICATE of REGISTRATION





CONTACT US:

(+91) 992-680-4526

(+91) 782-834-2091

harirajput0220@gmail.com

VISIT:

semssurgicalindustries.com

AUTHORIZED DISTRIBUTER:



ADDRESS:

❑ GOVINDPURA INDUSTRIAL AREA, BHOPAL, MADHYA PRADESH, INDIA

❑ 120, SEMS SURGICAL INDUSTRIES, AKVN ACHARPURA INDUSTRIAL AREA, BHOPAL, MADHYA PRADESH, INDIA



THE SYMBOL OF QUALITY & RELIABILITY

- Sems Surgical Industries Is A Manufacturer & Supplier Of Hospital Furniture & Equipments.
- Sems Surgical Industries Is A Licensed User Of The Trademark "Sems" & The Logo.
- Sems Surgical Industries Owns The Image Files In This Catalogue, Which Are Protected By Copyright And Other Restrictions.
- This Leaflet May Only Be Used For Offering, Purchasing, And Selling Our Products. No Part May Be Copied Or Reproduced In Any Way. In The Case Of Misuse, We Retain The Right To Recall Our Catalogues And Price Lists And To Take Legal Action.

- Sems Surgical Industries All Rights Reserved.